

TFCP2L1 Antibody


Catalog No: #25477

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TFCP2L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the amino terminus of human TFCP2L1.
Target Name	TFCP2L1
Other Names	Transcription factor CP2-like protein 1, LBP9, CRTR1
Accession No.	Swiss-Prot:Q9NZI6Gene ID:29842
Uniprot	Q9NZI6
GeneID	29842;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of TFCP2L1 in human colon tissue lysate with TFCP2L1 antibody at (A) 1 and (B) 2 ug/mL

Note: This product is for in vitro research use only