

ESRRB Antibody


Catalog No: #25479

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ESRRB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	19 amino acid peptide near the carboxy terminus of human ESRRB.
Target Name	ESRRB
Other Names	Estrogen-related receptor beta, ERR2, ERRb, ESRL2, NR3B2, DFNB35
Accession No.	Swiss-Prot:O95718Gene ID:2103
Uniprot	O95718
GeneID	2103;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of ESRRB in human heart tissue lysate with ESRRB antibody at 1 ug/mL in (A) the absence and (B) the presence of blocking peptide.

Note: This product is for in vitro research use only