

## NIBRIN Antibody


Catalog No: #25487

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NIBRIN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms
Specificity	Two alternatively spliced transcript isoforms of NIBRIN are known to exist.
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the carboxy terminus of human NIBRIN.
Target Name	NIBRIN
Other Names	Nijmegen breakage syndrome protein 1, NBN, AT-V1, AT-V2, ATV, NBS, NBS1, P95, Cell cycle regulatory protein p95
Accession No.	Swiss-Prot:O60934Gene ID:4683
Uniprot	O60934
GeneID	4683;
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Images



Note: This product is for in vitro research use only