

## RFX4 Antibody FITC Conjugated

Catalog No: #C01612F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RFX4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 125-175 735 derived from human RFX4
Conjugates	FITC
Target Name	RFX4
Other Names	NYD-SP10; Transcription factor RFX4; Regulatory factor X 4; Testis development protein NYD-SP10; RFX4
Accession No.	Swiss-Prot#Q33E94NCBI Gene ID5992
Uniprot	Q33E94
GeneID	5992;
Excitation Emission	494nm 518nm
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

May activate transcription by interacting directly with the X-box.

Note: This product is for in vitro research use only