

SOX21 Antibody FITC Conjugated

Catalog No: #C01632F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SOX21 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human SOX21
Conjugates	FITC
Target Name	SOX21
Other Names	Transcription factor SOX-21; SOX21; SOX-A; SOX25, SOXA
Accession No.	NCBI Gene ID11166
Uniprot	Q9Y651
GeneID	11166;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

May play a role as an activator of transcription of OPRM1.

Note: This product is for in vitro research use only