

## PHEX Antibody FITC Conjugated

Catalog No: #C01895F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PHEX Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 220-278 749 derived from human PHEX
Conjugates	FITC
Target Name	PHEX
Other Names	HYP; PEX; XLH; HPDR; HYP1; LXHR; HPDR1; Phosphate-regulating neutral endopeptidase; Metalloendopeptidase homolog PEX; Vitamin D-resistant hypophosphatemic rickets protein; X-linked hypophosphatemia protein; PHEX
Accession No.	Swiss-Prot#P78562NCBI Gene ID5251
Uniprot	P78562
GeneID	5251;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

Probably involved in bone and dentin mineralization and renal phosphate reabsorption.

Note: This product is for in vitro research use only