

ASCL2 Antibody FITC Conjugated

Catalog No: #C01928F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ASCL2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 193 derived from human ASCL2
Conjugates	FITC
Target Name	ASCL2
Other Names	ASH2; HASH2; MASH2; bHLHa45; Achaete-scute homolog 2; ASH-2; Class A basic helix-loop-helix protein 45; ASCL2
Accession No.	Swiss-Prot#Q99929NCBI Gene ID430
Uniprot	Q99929
GeneID	430;
Excitation Emission	494nm 518nm
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

AS-C proteins are involved in the determination of the neuronal precursors in the peripheral nervous system and the central nervous system.

Note: This product is for in vitro research use only