

human Fibrinogen Antibody HRP Conjugated

Catalog No: #C01982H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	human Fibrinogen Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC
Species Reactivity	Hu Ms Rt
Immunogen Description	Fibrinogen from human plasma
Conjugates	HRP
Target Name	human Fibrinogen
Other Names	Fib2; Fibrinogen alpha chain; FGA
Accession No.	Swiss-Prot#P02671NCBI Gene ID2243
Uniprot	P02671
GeneID	2243;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Cleaved by the protease thrombin to yield monomers which, together with fibrinogen beta (FGB) and fibrinogen gamma (FGG), polymerize to form an insoluble fibrin matrix. Fibrin has a major function in hemostasis as one of the primary components of blood clots. In addition, functions during the early stages of wound repair to stabilize the lesion and guide cell migration during re-epithelialization. Was originally thought to be essential for platelet aggregation, based on in vitro studies using anticoagulated blood. However, subsequent studies have shown that it is not absolutely required for thrombus formation in vivo. Enhances expression of SELP in activated platelets via an ITGB3-dependent pathway. Maternal fibrinogen is essential for successful pregnancy. Fibrin deposition is also associated with infection, where it protects against IFNG-mediated hemorrhage. May also facilitate the immune response via both innate and T-cell mediated pathways.

Note: This product is for in vitro research use only