

## CD28 Antibody FITC Conjugated

Catalog No: #C02234F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD28 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 20-45 218 derived from mouse CD28
Conjugates	FITC
Target Name	CD28
Other Names	T-cell-specific surface glycoprotein CD28; CD28
Accession No.	Swiss-Prot#P31041NCBI Gene ID12487
Uniprot	P31041
GeneID	12487;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.

Note: This product is for in vitro research use only