

CRK Monoclonal Antibody

Catalog No: #27127



Package Size: #27127-1 50ul #27127-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CRK Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3H7-E5-H8
Isotype	IgG2b
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Rb Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human CRK protein fragments expressed in E.coli.
Target Name	CRK
Other Names	p38, CRKII
Accession No.	Swiss-Prot:P05362Gene ID:3383
Uniprot	P05362
GeneID	3383;
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Application Details

Predicted MW: 34kd
Western blotting: 1:1000

Images



Western blot detection of CRK in LNCAP,A431

Background

The Crk-I and Crk-II forms differ in their biological activities. Crk-II has less transforming activity than Crk-I. Crk-II mediates attachment-induced

MAPK8 activation, membrane ruffling and cell motility in a Rac-dependent manner.

Note: This product is for in vitro research use only