

Syndecan 1 Antibody Biotin Conjugated

Catalog No: #C02324B

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syndecan 1 Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 275-310 310 derived from human CD138
Conjugates	Biotin
Target Name	CD138 Syndecan 1
Other Names	SDC; CD138; SYND1; syndecan; Syndecan-1; SDC1
Accession No.	Swiss-Prot#P18827NCBI Gene ID6382
Uniprot	P18827
GeneID	6382;
Excitation Emission	N A
Cell Localization	Cytoplasm, Membrane bound
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

Note: This product is for in vitro research use only