

FCGBP Antibody FITC Conjugated

Catalog No: #C02380F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FCGBP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 460-510 5405 derived from human FCGBP
Conjugates	FITC
Target Name	FCGBP
Other Names	FC(GAMMA)BP; IgGFc-binding protein; Fcgamma-binding protein antigen; FcgammaBP; FCGBP
Accession No.	Swiss-Prot#Q9Y6R7NCBI Gene ID8857
Uniprot	Q9Y6R7
GeneID	8857;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

May be involved in the maintenance of the mucosal structure as a gel-like component of the mucosa.

Note: This product is for in vitro research use only