

FGFR3 (Phospho-Tyr724) Antibody HRP Conjugated

Catalog No: #C02382H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | FGFR3 (Phospho-Tyr724) Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IHC-F ICC |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic phosphopeptide derived from human FGFR3 around the phosphorylation site of Tyr724 |
| Conjugates | HRP |
| Target Name | FGFR3 Tyr724 |
| Other Names | FGFR3 phospho Y724; FGFR3 phospho Tyr724; p-FGFR3 phospho Y724; ACH; CD 333; CD333; CD333 antigen; CEK 2; CEK2; FGFR 3; Fibroblast growth factor receptor 3 achondroplasia thanatophoric dwarfism; Fibroblast growth factor receptor 3; Heparin binding growth factor receptor; HSFGR3EX; Hydroxyaryl prote |
| Accession No. | NCBI Gene ID2261 |
| Uniprot | P22607 |
| GeneID | 2261; |
| Excitation Emission | N A |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC-F=1:50-200 ICC=1:50-200

Background

The encoded protein is synthesized mainly in corticotroph cells of the anterior pituitary where four cleavage sites are used; adrenocorticotrophin, essential for normal steroidogenesis and the maintenance of normal adrenal weight, and lipotropin beta are the major end products. In other tissues, including the hypothalamus, placenta, and epithelium, all cleavage sites may be used, giving rise to peptides with roles in pain and energy homeostasis, melanocyte stimulation, and immune modulation. These include several distinct melanotropins, lipotropins, and endorphins that are contained within the adrenocorticotrophin and beta-lipotropin peptides. Mutations in this gene have been associated with early onset obesity, adrenal insufficiency, and red hair pigmentation. Alternatively spliced transcript variants encoding the same protein have been described.

Note: This product is for in vitro research use only