

GIPR Antibody FITC Conjugated

Catalog No: #C02485F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GIPR Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 40-90 466 derived from human GIPR
Conjugates	FITC
Target Name	GIPR
Other Names	PGQTL2; Gastric inhibitory polypeptide receptor; GIP-R; Glucose-dependent insulinotropic polypeptide receptor; GIPR
Accession No.	Swiss-Prot#P48546NCBI Gene ID2696
Uniprot	P48546
GeneID	2696;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

This is a receptor for GIP. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only