

TIE1 Antibody FITC Conjugated

Catalog No: #C02532F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIE1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 850-900 1138 derived from human Tie1
Conjugates	FITC
Target Name	TIE1
Other Names	TIE; JTK14; Tyrosine-protein kinase receptor Tie-1; TIE1
Accession No.	Swiss-Prot#P35590NCBI Gene ID7075
Uniprot	P35590
GeneID	7075;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Transmembrane tyrosine-protein kinase that may modulate TEK TIE2 activity and contribute to the regulation of angiogenesis.

Note: This product is for in vitro research use only