

Rabbit anti-Goat IgG Secondary Antibody HRP conjugated

Catalog No: #L3042

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rabbit anti-Goat IgG Secondary Antibody HRP conjugated
Brief Description	Secondary Antibody
Clonality	Polyclonal
Applications	WB IHC ELISA
Species Reactivity	Goat
Immunogen Type	Peptide
Immunogen Description	The antiserum was developed in rabbit using the goat IgG as the immunogen.
Conjugates	HRP
Target Name	Goat IgG
Concentration	1mg/ml
Formulation	HRP-labeled Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 50% glycerol.
Storage	Store at -20 °C/1 year

Application Details

ELISA : 1:5000-1:20000
Western blotting: 1:5000-1:10000
Immunohistochemistry: 1:1000-5000

Note: This product is for in vitro research use only