

GPAT2 Antibody HRP Conjugated

Catalog No: #C02621H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPAT2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 60-110 795 derived from human GPAT2
Conjugates	HRP
Target Name	GPAT2
Other Names	CT123; Glycerol-3-phosphate acyltransferase 2, mitochondrial; GPAT-2; xGPAT1; GPAT2
Accession No.	Swiss-Prot#Q6NUI2NCBI Gene ID150763
Uniprot	Q6NUI2
GeneID	150763;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Esterifies acyl-group from acyl-ACP to the sn-1 position of glycerol-3-phosphate, an essential step in glycerolipid biosynthesis.

Note: This product is for in vitro research use only