

STAT6(Ab-645) Antibody

Catalog No: #21051



Package Size: #21051-1 50ul #21051-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	STAT6(Ab-645) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total STAT6 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.643~647 (P-A-T-I-K) derived from Human STAT6.
Target Name	STAT6
Other Names	IL-4-STAT; STAT6B; STAT6C
Accession No.	Swiss-Prot: P42226NCBI Protein: NP_003144.3
Uniprot	P42226
GeneID	6778;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 110kd
Western blotting: 1:500~1:1000
Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from HeLa cells using STAT6(Ab-645) Antibody #21051 and the same antibody preincubated with blocking peptide.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT6(Ab-645) Antibody #21051(left) or the same antibody preincubated with blocking peptide(right).

Background

Carries out a dual function: signal transduction and activation of transcription. Involved in interleukin-4 signalling.

Note: This product is for in vitro research use only