

Haemophilus influenza B Antibody HRP Conjugated

Catalog No: #C03354H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Haemophilus influenza B Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Haemophilus Influenzae
Crossing Reactivity	Haemophilus Influenzae
Immunogen Description	KLH conjugated synthetic peptide aa 100-150 339 derived from Haemophilus influenza B
Conjugates	HRP
Target Name	Haemophilus influenza B
Other Names	H. influenza B;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Haemophilus influenza B

Note: This product is for in vitro research use only