

HAPLN1 Antibody HRP Conjugated

Catalog No: #C03361H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HAPLN1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human HAPLN1
Conjugates	HRP
Target Name	HAPLN1
Other Names	Cartilage link protein; cartilage linking protein 1; CLP; CRTL1; Crtl1l; hyaluronan and proteoglycan link protein 1; Hyaluronan and proteoglycan link protein 1 precursor; LP; Proteoglycan link protein; HPLN1_HUMAN.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

HAPLN1

Note: This product is for in vitro research use only