

AQP5 Antibody FITC Conjugated

Catalog No: #C03486F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | AQP5 Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | Flow-Cyt IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 227-265 265 derived from mouse AQP5 |
| Conjugates | FITC |
| Target Name | Aquaporin 5 |
| Other Names | Aquaporin-5; AQP-5; Aqp5 |
| Accession No. | Swiss-Prot#Q9WTY4NCBI Gene ID11830 |
| Uniprot | Q9WTY4 |
| GeneID | 11830; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Cytoplasm |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Forms a water-specific channel. Implicated in the generation of saliva, tears, and pulmonary secretions.

Note: This product is for in vitro research use only