

Collagen VII Antibody FITC Conjugated

Catalog No: #C03520F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen VII Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 2944 derived from human Collagen VII
Conjugates	FITC
Target Name	Collagen VII
Other Names	EBD1; EBR1; EBDCT; Collagen alpha-1(VII) chain; Long-chain collagen; LC collagen; COL7A1
Accession No.	Swiss-Prot#Q02388NCBI Gene ID1294
Uniprot	Q02388
GeneID	1294;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Stratified squamous epithelial basement membrane protein that forms anchoring fibrils which may contribute to epithelial basement membrane organization and adherence by interacting with extracellular matrix (ECM) proteins such as type IV collagen.

Note: This product is for in vitro research use only