

GLP-1R Antibody PE Conjugated

Catalog No: #C03524P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLP-1R Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from rat GLP-1R:B 101-200/463
Conjugates	PE
Target Name	GLP-1R
Other Names	GLP-1R; GLP1Rc; Glucagon-like peptide 1 receptor; GLP-1 receptor; GLP-1-R; Glp1r
Accession No.	Swiss-Prot#P32301 NCBI Gene ID25051
Uniprot	P32301
GeneID	25051
Excitation Emission	480,565nm 578nm
Cell Localization	cytomembraneO cytoplasmO Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only