

GLP-1R Antibody HRP Conjugated

Catalog No: #C03524H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLP-1R Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 115-145 463 derived from rat GLP-1R
Conjugates	HRP
Target Name	GLP-1R
Other Names	GLP-1R; GLP1Rc; Glucagon-like peptide 1 receptor; GLP-1 receptor; GLP-1-R; Glp1r
Accession No.	Swiss-Prot#O35659NCBI Gene ID14652
Uniprot	O35659
GeneID	14652;
Excitation Emission	N A
Cell Localization	Extracellular, Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only