

## RAMP1 Antibody FITC Conjugated

Catalog No: #C03530F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RAMP1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 139-148 148 derived from mouse RAMP-1
Conjugates	FITC
Target Name	RAMP1
Other Names	9130218E19Rik; Receptor activity-modifying protein 1; Ramp1
Accession No.	Swiss-Prot#Q9WTJ5NCBI Gene ID51801
Uniprot	Q9WTJ5
GeneID	51801;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Transports the calcitonin gene-related peptide type 1 receptor (CALCRL) to the plasma membrane. Acts as a receptor for calcitonin-gene-related peptide (CGRP) together with CALCRL.

Note: This product is for in vitro research use only