

RAMP1 Antibody HRP Conjugated

Catalog No: #C03530H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RAMP1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 139-148 148 derived from mouse RAMP-1
Conjugates	HRP
Target Name	RAMP1
Other Names	9130218E19Rik; Receptor activity-modifying protein 1; Ramp1
Accession No.	Swiss-Prot#Q9WTJ5NCBI Gene ID51801
Uniprot	Q9WTJ5
GeneID	51801;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Transports the calcitonin gene-related peptide type 1 receptor (CALCRL) to the plasma membrane. Acts as a receptor for calcitonin-gene-related peptide (CGRP) together with CALCRL.

Note: This product is for in vitro research use only