

CK7 Antibody FITC Conjugated

Catalog No: #C03559F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CK7 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 265-315 496 derived from human of human CK7
Conjugates	FITC
Target Name	CK7
Other Names	K7; CK7; SCL; K2C7; Keratin, type II cytoskeletal 7; Cytokeratin-7; CK-7; Keratin-7; Sarcolectin; Type-II keratin Kb7; KRT7
Accession No.	Swiss-Prot#P08729NCBI Gene ID3855
Uniprot	P08729
GeneID	3855;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Blocks interferon-dependent interphase and stimulates DNA synthesis in cells. Involved in the translational regulation of the human papillomavirus type 16 E7 mRNA (HPV16 E7).

Note: This product is for in vitro research use only