

IDO2 Antibody HRP Conjugated

Catalog No: #C03598H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IDO2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F ICC
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 290-340 407 derived from human IDO2
Conjugates	HRP
Target Name	IDO2
Other Names	INDOL1; Indoleamine 2,3-dioxygenase 2; IDO-2; Indoleamine 2,3-dioxygenase-like protein 1; Indoleamine-pyrrole 2,3-dioxygenase-like protein 1; IDO2
Accession No.	Swiss-Prot#Q6ZQW0NCBI Gene ID169355
Uniprot	Q6ZQW0
GeneID	169355;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Catalyzes the first and rate-limiting step in the kynurenine pathway of tryptophan catabolism.

Note: This product is for in vitro research use only