

DRD3 Antibody FITC Conjugated

Catalog No: #C03671F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DRD3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 398-434 446 derived from mouse DRD3
Conjugates	FITC
Target Name	DRD3
Other Names	D3R; D(3) dopamine receptor; Dopamine D3 receptor; Drd3
Accession No.	Swiss-Prot#P30728NCBI Gene ID13490
Uniprot	P30728
GeneID	13490;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Dopamine receptor whose activity is mediated by G proteins which inhibit adenylyl cyclase. Promotes cell proliferation (By similarity).

Note: This product is for in vitro research use only