

EGFL7 Antibody HRP Conjugated

Catalog No: #C03677H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGFL7 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-80 275 derived from mouse EGFL7
Conjugates	HRP
Target Name	EGFL7
Other Names	Zneu1; VE-statin; Epidermal growth factor-like protein 7; EGF-like protein 7; Multiple epidermal growth factor-like domains protein 7; Multiple EGF-like domains protein 7; NOTCH4-like protein; Vascular endothelial statin; Egfl7; Megf7
Accession No.	Swiss-Prot#Q9QXT5NCBI Gene ID353156
Uniprot	Q9QXT5
GeneID	353156;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Regulates vascular tubulogenesis in vivo. Inhibits platelet-derived growth factor (PDGF)-BB-induced smooth muscle cell migration and promotes endothelial cell adhesion to the extracellular matrix and angiogenesis.

Note: This product is for in vitro research use only