

IL-2RB Antibody HRP Conjugated

Catalog No: #C03829H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-2RB Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 501-551 551 derived from human IL2 Receptor beta
Conjugates	HRP
Target Name	IL2 Receptor beta
Other Names	CD122; IL15RB; P70-75; Interleukin-2 receptor subunit beta; IL-2 receptor subunit beta; IL-2R subunit beta; IL-2RB; High affinity IL-2 receptor subunit beta; p75; IL2RB
Accession No.	Swiss-Prot#P14784NCBI Gene ID3560
Uniprot	P14784
GeneID	3560;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Receptor for interleukin-2. This beta subunit is involved in receptor mediated endocytosis and transduces the mitogenic signals of IL2.

Note: This product is for in vitro research use only