

CD162 Antibody FITC Conjugated

Catalog No: #C04029F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD162 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 260-300 412 derived from human PSGL-1 CD162
Conjugates	FITC
Target Name	CD162
Other Names	CLA; CD162; PSGL1; PSGL-1; P-selectin glycoprotein ligand 1; Selectin P ligand; SELPLG
Accession No.	Swiss-Prot#Q14242NCBI Gene ID6404
Uniprot	Q14242
GeneID	6404;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Background

A SLe(x)-type proteoglycan, which through high affinity, calcium-dependent interactions with E-, P- and L-selectins, mediates rapid rolling of leukocytes over vascular surfaces during the initial steps in inflammation. Critical for the initial leukocyte capture.

Note: This product is for in vitro research use only