

## CD30 Antibody FITC Conjugated

Catalog No: #C04143F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD30 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Pig Rab
Immunogen Description	KLH conjugated synthetic peptide aa 150-190 595 derived from human CD30
Conjugates	FITC
Target Name	CD30
Other Names	CD3; Ki-1; D1S166E; Tumor necrosis factor receptor superfamily member 8; CD3L receptor; Ki-1 antigen; Lymphocyte activation antigen CD3; TNFRSF8
Accession No.	Swiss-Prot#P28908NCBI Gene ID943
Uniprot	P28908
GeneID	943;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test IF=1:50-200

## Background

Receptor for TNFSF8 CD3L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.

Note: This product is for in vitro research use only