

## TFPI LACI Antibody FITC Conjugated

Catalog No: #C04172F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TFPI LACI Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 130-180 304 derived from human TFPI
Conjugates	FITC
Target Name	TFPI LACI
Other Names	EPI; TFI; LACI; TFPI1; Tissue factor pathway inhibitor; TFPI; Extrinsic pathway inhibitor; Lipoprotein-associated coagulation inhibitor
Accession No.	Swiss-Prot#P10646NCBI Gene ID7035
Uniprot	P10646
GeneID	7035;
Excitation Emission	494nm 518nm
Cell Localization	Secreted, Membrane bound
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Inhibits factor X (X(a)) directly and, in a Xa-dependent way, inhibits VIIa tissue factor activity, presumably by forming a quaternary Xa LACI VIIa TF complex. It possesses an antithrombotic action and also the ability to associate with lipoproteins in plasma.

Note: This product is for in vitro research use only