

## CXCL5 Antibody HRP Conjugated

Catalog No: #C04181H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CXCL5 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-110 132 derived from mouse CXCL5
Conjugates	HRP
Target Name	CXCL5
Other Names	LIX; GCP-2; Scyb5; Scyb6; ENA-78; AMCF-II; C-X-C motif chemokine 5; Cytokine LIX; Small-inducible cytokine B5; Cxcl5
Accession No.	Swiss-Prot#P50228NCBI Gene ID20311
Uniprot	P50228
GeneID	20311;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

May participate in the recruitment of inflammatory cells by injured or infected tissue. GCP-2(1-78) and, more potent, GCP-2(9-78) attract neutrophils and are involved in neutrophil activation.

Note: This product is for in vitro research use only