

LAMP3 Antibody FITC Conjugated

Catalog No: #C04187F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LAMP3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 320-370 416 derived from human LAMP-3
Conjugates	FITC
Target Name	LAMP3
Other Names	LAMP; CD208; DCLAMP; LAMP-3; TSC403; DC LAMP; DC-LAMP; Lysosome-associated membrane glycoprotein 3; Lysosomal-associated membrane protein 3; DC-lysosome-associated membrane glycoprotein; Protein TSC403; LAMP3
Accession No.	Swiss-Prot#Q9UQV4NCBI Gene ID27074
Uniprot	Q9UQV4
GeneID	27074;
Excitation Emission	494nm 518nm
Cell Localization	Lumenal
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1ug/Test IF=1:50-200

Background

May play a role in dendritic cell function and in adaptive immunity.

Note: This product is for in vitro research use only