

GCSF Receptor Antibody FITC Conjugated

Catalog No: #C04197F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GCSF Receptor Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HuB MsB Rt
Immunogen Description	KLH conjugated synthetic peptide aa 125-175 863 derived from human G-CSFR
Conjugates	FITC
Target Name	G-CSFR
Other Names	CD114; GCSFR; Granulocyte colony-stimulating factor receptor; G-CSF receptor; G-CSF-R; CSF3R
Accession No.	Swiss-Prot#Q99062NCBI Gene ID1441
Uniprot	Q99062
GeneID	1441;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200B

Background

Receptor for granulocyte colony-stimulating factor (CSF3), essential for granulocytic maturation. Plays a crucial role in the proliferation, differentiation and survival of cells along the neutrophilic lineage. In addition it may function in some adhesion or recognition events at the cell surface.

Note: This product is for in vitro research use only