

## C3a Receptor C3aR Antibody FITC Conjugated

Catalog No: #C04389F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C3a Receptor C3aR Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HuB MsB RtB B B
Immunogen Description	KLH conjugated synthetic peptide aa 439-482 482 derived from human C3a Receptor C3aR
Conjugates	FITC
Target Name	C3a Receptor C3aR
Other Names	AZ3B; C3AR; HNFAG09; C3a anaphylatoxin chemotactic receptor; C3a-R; C3AR1; C3R1
Accession No.	Swiss-Prot#Q16581NCBI Gene ID719
Uniprot	Q16581
GeneID	719;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200B

## Background

Receptor for the chemotactic and inflammatory peptide anaphylatoxin C3a. This receptor stimulates chemotaxis, granule enzyme release and superoxide anion production.

Note: This product is for in vitro research use only