

RPS6 (Phospho-Ser240 + Ser244) Antibody FITC Conjugated

Catalog No: #C04670F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RPS6 (Phospho-Ser240 + Ser244) Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HuB MsB RtB
Immunogen Description	KLH conjugated synthetic phosphopeptide aa 210-249 249 derived from human S6 Ribosomal Protein around the phosphorylation site of Ser240 244
Conjugates	FITC
Target Name	RPS6 Ser240 + Ser244
Other Names	S6; 40S ribosomal protein S6; Phosphoprotein NP33; RPS6; OK SW-cl.2
Accession No.	Swiss-Prot#P62753NCBI Gene ID6194
Uniprot	P62753
GeneID	6194;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200B

Background

May play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA.

Note: This product is for in vitro research use only