

## IL18BP Antibody FITC Conjugated

Catalog No: #C05106F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL18BP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 194 derived from human IL18BP
Conjugates	FITC
Target Name	IL18BP
Other Names	IL18BPα; Interleukin-18-binding protein; IL-18BP; Tadekinig-alfa; IL18BP
Accession No.	Swiss-Prot#O95998NCBI Gene ID10068
Uniprot	O95998
GeneID	10068;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Note: This product is for in vitro research use only