

tubulin Beta Antibody HRP Conjugated

Catalog No: #C05295H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	tubulin Beta Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 444 derived from human tubulin Beta
Conjugates	HRP
Target Name	tubulin Beta
Other Names	M4; TUBB1; TUBB5; CDCBM6; OK SW-cl.56; Tubulin beta chain; Tubulin beta-5 chain; TUBB
Accession No.	Swiss-Prot#P07437NCBI Gene ID203068
Uniprot	P07437
GeneID	203068;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-1:50000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Note: This product is for in vitro research use only