

CD1d1 Antibody PE Conjugated

Catalog No: #C05458P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD1d1 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from mouse CD1d1
Conjugates	PE
Target Name	CD1d1
Other Names	Cd1a; Cd1d; CD1.1; Ly-38; AI747460; Antigen-presenting glycoprotein CD1d1; CD1d.1; Cd1d1
Accession No.	Swiss-Prot#P11609NCBI Gene ID12479
Uniprot	P11609
GeneID	12479;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Antigen-presenting protein that binds self and non-self glycolipids and presents them to T-cell receptors on natural killer T-cells.

Note: This product is for in vitro research use only