

Mouse IgG1 Isotype Control, PE Conjugated mAb

Catalog No: #28268

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Mouse IgG1 Isotype Control, PE Conjugated mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A1B1
Isotype	IgG1
Applications	FC
Conjugates	PE
Target Name	IgG1
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Vol.per.Test: 10 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically add 10ul of this reagent to 100 μ l of experimental sample with 1 X 10⁶ cells per test. Please refer to the detailed protocol when you perform a test.

Images



Human peripheral blood lymphocytes analyzed with PE
Mouse IgG1 Isotype Control mAb

Background

This antibody can be used as a mouse IgG1 isotype control in flow cytometry and other applications.

Note: This product is for in vitro research use only