

Mouse IgG2b Isotype Control, PE Conjugated mAb

Catalog No: #28280

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


| | |
|--------------|---|
| Product Name | Mouse IgG2b Isotype Control, PE Conjugated mAb |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 4A2B5 |
| Isotype | IgG2b |
| Applications | FC |
| Conjugates | PE |
| Target Name | IgG2b |
| Formulation | Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Vol.per.Test: 10 &mu;l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically add 10ul of this reagent to 100μl of experimental sample with 1 X 10⁶ cells per test. Please refer to the detailed protocol when you perform a test.

Images



Human peripheral blood lymphocytes analyzed with PE
Mouse IgG2b Isotype Control mAb

Background

This antibody can be used as a mouse IgG2b isotype control in flow cytometry and other applications.

Note: This product is for in vitro research use only