

LHR Antibody HRP Conjugated

Catalog No: #C06407H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LHR Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 210-260 699 derived from human LHR CGR
Conjugates	HRP
Target Name	LHR
Other Names	HHG; LHR; LCGR; LGR2; ULG5; LHRHR; LSH-R; LH CGR; LH CG-R; Lutropin-choriogonadotropic hormone receptor; Luteinizing hormone receptor; LHCGR
Accession No.	Swiss-Prot#P22888NCBI Gene ID3973
Uniprot	P22888
GeneID	3973;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Receptor for lutropin-choriogonadotropic hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

Note: This product is for in vitro research use only