

## Raf1(Ab-338) Antibody

Catalog No: #21202

Package Size: #21202-1 50ul #21202-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Raf1(Ab-338) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total Raf1 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.336~340 (R-D-S-S-Y) derived from Human RAF.
Target Name	Raf1
Other Names	C-RAF; C-Raf; CRAF
Accession No.	Swiss-Prot: P04049 NCBI Protein: NP_002871.1
Uniprot	P04049
GeneID	5894;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 73kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100

## Images



Western blot analysis of extracts from 293 cells using Raf1(Ab-338) Antibody #21202 and the same antibody preincubated with blocking peptide.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Raf1(Ab-338) Antibody #21202(left) or the same antibody preincubated with blocking peptide(right).

## Background

Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. Part of the Ras-dependent signaling pathway from receptors to the nucleus. Protects cells from apoptosis mediated by STK3.

Adams DG,et al. (2005)J Biol Chem. 280(52):42644-54.

Miura S, et al. (2003) Biochem Biophys Res Commun. 306(4):924-9.

Noser JA,et al.(2007)Mol Ther.15(8):1531-6.

Note: This product is for in vitro research use only