

RNASE13 Antibody

Catalog No: #32146

Package Size: #32146-1 50ul #32146-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RNASE13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total RNASE13 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human RNASE13 (NP_001012264.1).
Target Name	RNASE13
Other Names	RNASE13;HEL-S-86p;RAL1
Accession No.	Uniprot:Q5GAN3GeneID:440163
Uniprot	Q5GAN3
GeneID	440163
SDS-PAGE MW	18kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:100 - 1:200


Images




Immunohistochemistry of paraffin-embedded human prostate using RNASE13 antibody.



Western blot analysis of extracts of Raji cells, using RNASE13 antibody.



Immunohistochemistry of paraffin-embedded mouse lung using RNASE13 antibody.



Immunohistochemistry of paraffin-embedded human stomach using RNASE13 antibody.

Background

Does not exhibit any ribonuclease activity.

Note: This product is for in vitro research use only