

RNF19B Antibody HRP Conjugated

Catalog No: #C07840H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNF19B Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ELISA IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human RNF19B
Conjugates	HRP
Target Name	RNF19B
Other Names	E3 ubiquitin-protein ligase RNF19B; IBR domain containing 3; IBR domain-containing protein 3; IBRDC3; Natural killer lytic-associated molecule; NKLAM; RING finger protein 19B; RN19B_HUMAN; rnf19b.
Accession No.	NCBI Gene ID127544
Uniprot	Q6ZMZ0
GeneID	127544;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P□1:100-500

IHC-F□1:100-500

ELISA□1:5000-10000

Background

E3 ubiquitin-protein ligase which accepts ubiquitin from E2 ubiquitin-conjugating enzymes UBE2L3 and UBE2L6 in the form of a thioester and then directly transfers the ubiquitin to targeted substrates, such as UCKL1. Involved in the cytolytic activity of natural killer cells and cytotoxic T-cells.

Note: This product is for in vitro research use only