

LETMD1 Antibody

Catalog No: #32624



Package Size: #32624-1 50ul #32624-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LETMD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total LETMD1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human LETMD1.
Target Name	LETMD1
Other Names	HCCR1; HCCR-1; HCCR-2; HCRR-2; 1110019O13Rik
Accession No.	Swiss-Prot:Q6P1Q0NCBI Gene ID:25875
Uniprot	Q6P1Q0
GeneID	25875;
SDS-PAGE MW	42KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using LETMD1 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded human liver damage using LETMD1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat kidney using LETMD1 at dilution of 1:100 (40x lens).

Background

This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene.

Note: This product is for in vitro research use only