

MSN Antibody

Catalog No: #32645



Package Size: #32645-1 50ul #32645-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	MSN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total MSN protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human Moesin (NP_002435.1).
Target Name	MSN
Other Names	MSN;HEL70;IMD50;moesin
Accession No.	Uniprot:P26038GeneID:4478
Uniprot	P26038
GeneID	4478
SDS-PAGE MW	88kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using Moesin antibody.

## Background

---

Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement.

---

Note: This product is for in vitro research use only